

AKAP12 antibody

Product Information

Catalog No.:	FNab00252
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Anchoring protein that mediates the subcellular compartmentation of protein kinase A(PKA) and protein kinase C(PKC).

Immunogen information


Immunogen:	A kinase(PRKA) anchor protein 12
Synonyms:	A-kinase anchor protein 12 (AKAP-12) A-kinase anchor protein 250 kDa (AKAP 250) Gravin Myasthenia gravis autoantigen AKAP12 AKAP250
Observed MW:	191 kDa
Uniprot ID :	Q02952

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-200
Image:	



Immunohistochemistry of paraffin-embedded human colon tissue slide using FNab00252(AKAP12 Antibody) at dilution of 1:50



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab00252(AKAP12 Antibody) at dilution of 1:600