

## AKAP3 antibody

### Product Information

Catalog No.:	FNab00254
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May function as a regulator of both motility-and head-associated functions such as capacitation and the acrosome reaction.

### Immunogen information


Immunogen:	A kinase(PRKA) anchor protein 3
Synonyms:	A-kinase anchor protein 3 (AKAP-3) A-kinase anchor protein 110 kDa (AKAP 110) Cancer/testis antigen 82 (CT82) Fibrous sheath protein of 95 kDa (FSP95) Fibrousheathin I Fibrousheathin-1 Protein kinase A-anchoring protein 3 (PRKA3) Sperm oocyte-binding protein AKAP3 AKAP110 SOB1
Observed MW:	95 kDa
Uniprot ID :	O75969

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:20-200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab00254(AKAP3 antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00254(AKAP3 antibody) at dilution of 1:5000