

AKAP4 antibody

Product Information

Catalog No.:	FNab00256
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Major structural component of sperm fibrous sheath. Plays a role in sperm motility.

Immunogen information


Immunogen:	A kinase(PRKA) anchor protein 4
Synonyms:	A-kinase anchor protein 4 (AKAP-4) A-kinase anchor protein 82 kDa (AKAP 82, hAKAP82) Major sperm fibrous sheath protein (HI) Protein kinase A-anchoring protein 4 (PRKA4) AKAP4 AKAP82
Observed MW:	94 kDa
Uniprot ID :	Q5JQC9

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab00256(AKAP4 Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab00256(AKAP4 Antibody) at dilution of 1:1500