

ALLC antibody

Product Information

Catalog No.:	FNab00316
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The function of this enzyme is unclear as allantoicase activity is not known to exist in mammals.

Immunogen information


Immunogen:	allantoicase
Synonyms:	Probable inactive allantoicase Allantoate amidinohydrolase ALLC
Observed MW:	46 kDa
Uniprot ID :	Q8N6M5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00316(ALLC Antibody) at dilution of 1:50



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab00316(ALLC antibody) at dilution of 1:500