

anti- Alpha actin antibody

Product Information

Catalog No.:	FNab00320
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	actin, alpha, cardiac muscle 1
Synonyms:	ACTC, ACTC1, ACTC1 cardiac muscle 1, Actin, alpha cardiac muscle 1, actin, alpha, cardiac muscle 1, Alpha actin, Alpha cardiac actin, CMD1R, CMH11
Observed MW:	42 kDa
UniprotID :	P68032

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB,IP
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-alpha Actin (IP:FNab00320, 4ug;
Detection:FNab00320 1:800) with mouse heart
tissue lysate 3200ug.



mouse heart tissue were subjected to SDS PAGE
followed by western blot with FNab00320(ACTC1
antibody) at dilution of 1:800