

## CTNNA2 antibody

### Product Information

Catalog No.:	FNab00331
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	catenin(cadherin-associated protein), alpha 2
Synonyms:	Catenin alpha-2 Alpha N-catenin Alpha-catenin-related protein CTNNA2 CAPR
Observed MW:	102 kDa
Uniprot ID :	P26232

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:1000-1:10000; IP: 1:200-1:2000; IF: 1:20-1:200
Image:	



IP Result of anti-alpha N Catenin (IP:FNab00331, 4ug; Detection:FNab00331 1:500) with SH-SY5Y cells lysate 1600ug.



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab00331(CTNNA2 antibody) at dilution of 1:1200