

## SGCA antibody

### Product Information

Catalog No.:	FNab00342
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	sarcoglycan, alpha(50kDa dystrophin-associated glycoprotein)
Synonyms:	Alpha-sarcoglycan (Alpha-SG) 50 kDa dystrophin-associated glycoprotein (50DAG) Adhalin Dystroglycan-2 SGCA ADL DAG2
Observed MW:	43 kDa
Uniprot ID :	Q16586

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab00342(Alpha-sarcoglycan antibody) at dilution of 1:1000