

ALYREF antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab00354 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information

| | |
|--------------|---|
| Immunogen: | THO complex 4 |
| Synonyms: | THO complex subunit 4 (Tho4) Ally of AML-1 and LEF-1 Aly/REF export factor Transcriptional coactivator Aly/REF bZIP-enhancing factor BEF ALYREF ALY BEF THOC4 |
| Observed MW: | 32 kDa |
| Uniprot ID : | Q86V81 |

Application

| | |
|--------------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:200-1:1000 |
| Image: | |



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00354(ALY Antibody) at dilution of 1:300