

PRKAG2 antibody

Product Information

Catalog No.:	FNab00380
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	protein kinase, AMP-activated, gamma 2 non-catalytic subunit
Synonyms:	5'-AMP-activated protein kinase subunit gamma-2 (AMPK gamma2, AMPK subunit gamma-2) H91620p PRKAG2
Observed MW:	63 kDa
Uniprot ID :	Q9UGJ0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab00380 (PRKAG2 Antibody) at dilution of 1:50



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab00380 (PRKAG2 Antibody) at dilution of 1:600