

## ANGPTL5 antibody

### Product Information

Catalog No.:	FNab00400
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	angiopoietin-like 5
Synonyms:	Angiopoietin-related protein 5 Angiopoietin-like protein 5 ANGPTL5
Observed MW:	44 kDa
Uniprot ID :	Q86XS5

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab00400(ANGPTL5 antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00400(ANGPTL5 antibody) at dilution of 1:500