
ANKMY2 antibody

Product Information

Catalog No.:	FNab00403
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in the trafficking of signaling proteins to the cilia.

Immunogen information

Immunogen:	ankyrin repeat and MYND domain containing 2
Synonyms:	Ankyrin repeat and MYND domain-containing protein 2 ANKMY2
Observed MW:	49 kDa
Uniprot ID :	Q8IV38

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00403(ANKMY2 antibody) at dilution of 1:500