

ANKRD11 antibody

Product Information

Catalog No.:	FNab00405
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


ANKRD11, also named as ANCO1, may recruit HDACs to the p160 coactivators/nuclear receptor complex to inhibit ligand-dependent transactivation. This antibody is specific to ANKRD11. ANKRD11 exists 150-170 kDa isoform.

Immunogen information

Immunogen:	ankyrin repeat domain 11
Synonyms:	Ankyrin repeat domain-containing protein 11 Ankyrin repeat-containing cofactor 1 ANKRD11 ANCO1
Observed MW:	160 kDa
Uniprot ID :	Q6UB99

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab00405(ANKRD11 antibody) at dilution of 1:500