

ANKRD13C antibody

Product Information

Catalog No.:	FNab00408
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a molecular chaperone for G protein-coupled receptors, regulating their biogenesis and exit from the ER.

Immunogen information


Immunogen:	ankyrin repeat domain 13C
Synonyms:	Ankyrin repeat domain-containing protein 13C ANKRD13C
Observed MW:	61 kDa
Uniprot ID :	Q8N6S4

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab00408(ANKRD13C Antibody) at dilution of 1:50



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with FNab00408(ANKRD13C Antibody) at dilution of 1:600