

ANKRD26 antibody

Product Information

Catalog No.:	FNab00411
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ankyrin repeat domain 26
Synonyms:	Ankyrin repeat domain-containing protein 26 ANKRD26 KIAA1074
Observed MW:	196 kDa
Uniprot ID :	Q9UPS8

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, IF, IP
Recommended dilution:	IP: 1:10-1:100; IHC: 1:20-1:200; IF: 1:200-1:2000
Image:	



Immunohistochemistry of paraffin-embedded human ovary using FNab00411(anti- ANKRD26 antibody) at dilution of 1:50



IP Result of anti-ANKRD26 (IP:FNab00411, 3ug; Detection:FNab00411 1:600) with HeLa cells lysate 4000ug.