

ANKRD7 antibody

Product Information

Catalog No.:	FNab00421
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	ankyrin repeat domain 7
Synonyms:	Ankyrin repeat domain-containing protein 7 Testis-specific protein TSA806 ANKRD7
Observed MW:	29 kDa
Uniprot ID :	Q92527

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab00421(ANKRD7 Antibody) at dilution of 1:300