

ANXA10 antibody

Product Information

Catalog No.:	FNab00427
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined.

Immunogen information

Immunogen:	annexin A10
Synonyms:	Annexin A10 Annexin-10 Annexin-14 ANXA10 ANX14
Observed MW:	37 kDa
Uniprot ID :	Q9UJ72

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



U-251MG cells were subjected to SDS PAGE followed by western blot with FNab00427(ANXA10 Antibody) at dilution of 1:2000