

## ANXA8 antibody

### Product Information

Catalog No.:	FNab00440
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	annexin A8-like 2
Synonyms:	Annexin A8 Annexin VIII Annexin-8 Vascular anticoagulant-beta (VAC-beta) ANXA8 ANX8
Observed MW:	37 kDa
Uniprot ID :	P13928

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human skin tissue slide using FNab00440(Annexin VIII antibody) at dilution of 1:200



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00440(ANXA8 antibody) at dilution of 1:100