

ANO8 antibody

Product Information

Catalog No.:	FNab00444
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Does not exhibit calcium-activated chloride channel(CaCC) activity.

Immunogen information

Immunogen:	anoctamin 8
Synonyms:	Anoctamin-8 Transmembrane protein 16H ANO8 KIAA1623 TMEM16H
Observed MW:	136 kDa
Uniprot ID :	Q9HCE9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00444(ANO8 antibody) at dilution of 1:500