

AP1G2 antibody

Product Information

Catalog No.:	FNab00458
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function in protein sorting in late endosomes or multivesicular bodies(MVBs). Involved in MVB-assisted maturation of hepatitis B virus(HBV).

Immunogen information

Immunogen:	adaptor-related protein complex 1, gamma 2 subunit
Synonyms:	AP-1 complex subunit gamma-like 2 Gamma2-adaptin (G2ad) AP1G2
Observed MW:	70 kDa, 87 kDa
Uniprot ID :	O75843

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



THP-1 cells were subjected to SDS PAGE followed by western blot with FNab00458(AP1G2 Antibody) at dilution of 1:500