

## APBA2 antibody

### Product Information

Catalog No.:	FNab00471
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Putative function in synaptic vesicle exocytosis by binding to STXBP1, an essential component of the synaptic vesicle exocytotic machinery. May modulate processing of the beta-amyloid precursor protein(APP) and hence formation of beta-APP.

### Immunogen information

Immunogen:	amyloid beta(A4) precursor protein-binding, family A, member 2
Synonyms:	Amyloid-beta A4 precursor protein-binding family A member 2 Adapter protein X11beta Neuron-specific X11L protein Neuronal Munc18-1-interacting protein 2 (Mint-2) APBA2 MINT2 X11L
Observed MW:	120 kDa
Uniprot ID :	Q99767

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00471(APBA2 antibody) at dilution of 1:500