

ANAPC1 antibody

Product Information

Catalog No.:	FNab00474
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	anaphase promoting complex subunit 1
Synonyms:	Anaphase-promoting complex subunit 1 (APC1) Cyclosome subunit 1 Mitotic checkpoint regulator Testis-specific gene 24 protein ANAPC1 TSG24
Observed MW:	200 kDa
Uniprot ID :	Q9H1A4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab00474(ANAPC1 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00474(ANAPC1 antibody) at dilution of 1:1500