

APLN antibody

Product Information

Catalog No.:	FNab00478
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	apelin
Synonyms:	Apelin APJ endogenous ligand Apelin-36 Apelin-31 Apelin-28 Apelin-13 APLN APEL
Observed MW:	None
Uniprot ID :	Q9ULZ1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse heart tissue using FNab00478(anti- Apelin antibody) at dilution of 1:200