

## APOBEC2 antibody

### Product Information

Catalog No.:	FNab00492
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probable C to U editing enzyme whose physiological substrate is not yet known. Does not display detectable apoB mRNA editing. Has a low intrinsic cytidine deaminase activity. May play a role in the epigenetic regulation of gene expression through the process of active DNA demethylation.

### Immunogen information


Immunogen:	apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 2
Synonyms:	C->U-editing enzyme APOBEC-2 mRNA(cytosine(6666)) deaminase 2 APOBEC2
Observed MW:	26 kDa
Uniprot ID :	Q9Y235

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab00492(APOBEC2 antibody) at dilution of 1:100



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab00492(APOBEC2 antibody) at dilution of 1:500