

APOD antibody

Product Information

Catalog No.:	FNab00499
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


APOD occurs in the macromolecular complex with lecithin-cholesterol acyltransferase. It is probably involved in the transport and binding of bilin. Appears to be able to transport a variety of ligands in a number of different contexts.

Immunogen information


Immunogen:	apolipoprotein D
Synonyms:	Apolipoprotein D (Apo-D, ApoD) APOD
Observed MW:	22-33 kDa
Uniprot ID :	P05090

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab00499(APOD antibody) at dilution of 1:100



human serum tissue were subjected to SDS PAGE followed by western blot with FNab00499(APOD antibody) at dilution of 1:1000