

## APOF antibody

### Product Information

Catalog No.:	FNab00501
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Minor apolipoprotein that associates with LDL. Inhibits cholesteryl ester transfer protein(CETP) activity and appears to be an important regulator of cholesterol transport. Also associates to a lesser degree with VLDL, Apo-AI and Apo-AII.

### Immunogen information


Immunogen:	apolipoprotein F
Synonyms:	Apolipoprotein F (Apo-F) Lipid transfer inhibitor protein (LTIP) APOF
Observed MW:	29 kDa
Uniprot ID :	Q13790

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNa00501(APOF antibody) at dilution of 1:50



human blood tissue were subjected to SDS PAGE followed by western blot with FNa00501(APOF antibody) at dilution of 1:1000