

## APOL1 antibody

### Product Information

Catalog No.:	FNab00503
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5E2
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in lipid exchange and transport throughout the body. May participate in reverse cholesterol transport from peripheral cells to the liver.

### Immunogen information


Immunogen:	apolipoprotein L, 1
Synonyms:	Apolipoprotein L1 Apolipoprotein L (Apo-L, ApoL) Apolipoprotein L-I (ApoL-I) APOL1 APOL
Observed MW:	44 kDa
Uniprot ID :	O14791

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:10000-1:50000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab00503(APOL1 antibody) at dilution of 1:50



1 ug human plasma tissue was subjected to SDS PAGE followed by western blot with FNab00503(APOL1 Antibody) at dilution of 1:50000