

## APOA2 antibody

### Product Information

Catalog No.:	FNab00507
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


The predicted MW of ApoA2 is 11 kDa, while the mature form is smaller (7-10 kDa) since the signal peptide was cleaved. In humans, most of circulating apoA2 exist as dimer.

### Immunogen information

Immunogen:	apolipoprotein A-II
Synonyms:	Apolipoprotein A-II (Apo-AII, ApoA-II) Apolipoprotein A2 Proapolipoprotein A-II (ProapoA-II) Truncated apolipoprotein A-II Alternative names: Apolipoprotein A-II(1-76) APOA2
Observed MW:	15 kDa, 11 kDa
Uniprot ID :	P02652

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis tissue slide using FNab00507(APOA2 Antibody) at dilution of 1:50



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab00507(APOA2 antibody) at dilution of 1:1000