

## NAE1 antibody

### Product Information

Catalog No.:	FNab00514
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	NEDD8 activating enzyme E1 subunit 1
Synonyms:	NEDD8-activating enzyme E1 regulatory subunit Amyloid beta precursor protein-binding protein 1, 59 kDa (APP-BP1) Amyloid protein-binding protein 1 Proto-oncogene protein 1 NAE1 APPBP1
Observed MW:	60 kDa
Uniprot ID :	Q13564

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-APPBP1 (IP:FNab00514, 4ug;  
Detection:FNab00514 1:500) with mouse kidney  
tissue lysate 4000ug.