

ARD1A antibody

Product Information

Catalog No.:	FNab00526
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ARD1 homolog A, N-acetyltransferase
Synonyms:	N-terminal amino-acid N(alpha)-acetyltransferase NatA ARD1A
Observed MW:	26 kDa
Uniprot ID :	Q6P4J0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab00526(ARD1A antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab00526(ARD1A antibody) at dilution of 1:500