

ARFIP1 antibody

Product Information

Catalog No.:	FNab00539
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Putative target protein of ADP-ribosylation factor.

Immunogen information

Immunogen:	ADP-ribosylation factor interacting protein 1
Synonyms:	Arfaptin-1 ADP-ribosylation factor-interacting protein 1 ARFIP1
Observed MW:	44 kDa
Uniprot ID :	P53367

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas using FNab00539(ARFIP1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00539(ARFIP1 antibody) at dilution of 1:500