

ARFRP1 antibody

Product Information

Catalog No.:	FNab00541
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Trans-Golgi-associated GTPase that regulates protein sorting. Controls the targeting of ARL1 and its effector to the trans-Golgi. Required for the lipidation of chylomicrons in the intestine and required for VLDL lipidation in the liver.

Immunogen information

Immunogen:	ADP-ribosylation factor related protein 1
Synonyms:	ADP-ribosylation factor-related protein 1 (ARF-related protein 1, ARP) ARFRP1 ARP1
Observed MW:	22-25 kDa
Uniprot ID :	Q13795

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



293T cells were subjected to SDS PAGE followed by western blot with FNab00541 (ARFRP1 antibody) at dilution of 1:500