

ARHGAP15 antibody

Product Information

Catalog No.:	FNab00547
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state. Has activity toward RAC1. Overexpression results in an increase in actin stress fibers and cell contraction.

Immunogen information


Immunogen:	Rho GTPase activating protein 15
Synonyms:	Rho GTPase-activating protein 15 ArhGAP15 Rho-type GTPase-activating protein 15 ARHGAP15
Observed MW:	65-70 kDa
Uniprot ID :	Q53QZ3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab00547(ARHGAP15 antibody) at dilution of 1:50



human liver tissue were subjected to SDS PAGE followed by western blot with FNab00547(ARHGAP15 antibody) at dilution of 1:500