

ARHGAP21 antibody

Product Information

Catalog No.:	FNab00550
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Functions as a GTPase-activating protein(GAP) for RHOA and CDC42. Downstream partner of ARF1 which may control Golgi apparatus structure and function. Also required for CTNNA1 recruitment to adherens junctions.

Immunogen information


Immunogen:	Rho GTPase activating protein 21
Synonyms:	Rho GTPase-activating protein 21 Rho GTPase-activating protein 10 Rho-type GTPase-activating protein 21 ARHGAP21 ARHGAP10 KIAA1424
Observed MW:	250 kDa
Uniprot ID :	Q5T5U3

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab00550 (ARHGAP21 Antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab00550 (ARHGAP21 Antibody) at dilution of 1:600