

## ARHGEF10 antibody

### Product Information

Catalog No.:	FNab00559
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May play a role in developmental myelination of peripheral nerves.

### Immunogen information

Immunogen:	Rho guanine nucleotide exchange factor(GEF) 10
Synonyms:	Rho guanine nucleotide exchange factor 10 ARHGEF10 KIAA0294
Observed MW:	68 kDa, 75 kDa, 148-152 kDa
Uniprot ID :	O15013

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab00559 (ARHGEF10 Antibody ) at dilution of 1:200