

## NUAK1 antibody

### Product Information

Catalog No.:	FNab00570
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	NUAK family, SNF1-like kinase, 1
Synonyms:	NUAK family SNF1-like kinase 1 AMPK-related protein kinase 5 (ARK5) Omphalocele kinase 1 NUAK1 ARK5 KIAA0537 OMPHK1
Observed MW:	74 kDa
Uniprot ID :	O60285

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab00570 (NUAK1 Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00570 (NUAK1 Antibody) at dilution of 1:600