

ARL1 antibody

Product Information

Catalog No.:	FNab00571
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


GTP-binding protein that has very low efficiency as allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Can activate phospholipase D with very low efficiency. Important for normal function of the Golgi apparatus.

Immunogen information

Immunogen:	ADP-ribosylation factor-like 1
Synonyms:	ADP-ribosylation factor-like protein 1 ARL1
Observed MW:	20 kDa
Uniprot ID :	P40616

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab00571(ARL1 Antibody) at dilution of 1:600