

ARL8A antibody

Product Information

Catalog No.:	FNab00581
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in lysosomes motility. Alternatively, may play a role in chromosome segregation(By similarity).

Immunogen information


Immunogen:	ADP-ribosylation factor-like 8A
Synonyms:	ADP-ribosylation factor-like protein 8A ADP-ribosylation factor-like protein 10B Novel small G protein indispensable for equal chromosome segregation 2 ARL8A ARL10B GIE2
Observed MW:	22 kDa
Uniprot ID :	Q96BM9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000
Image:	



IP Result of anti-ARL8A (IP:FNab00581, 4ug;
Detection:FNab00581 1:600) with mouse brain
tissue lysate 3440ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with
FNab00581(ARL8A Antibody) at dilution of
1:600