

ARMC8 antibody

Product Information

Catalog No.:	FNab00587
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	armadillo repeat containing 8
Synonyms:	Armadillo repeat-containing protein 8 ARMC8
Observed MW:	75.5 kDa
Uniprot ID :	Q8IUR7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:10000; IHC: 1:50-1:500
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab00587(ARMC8 Antibody) at dilution of 1:6000