

ARMCX1 antibody

Product Information

Catalog No.:	FNab00588
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	armadillo repeat containing, X-linked 1
Synonyms:	Armadillo repeat-containing X-linked protein 1 ARM protein lost in epithelial cancers on chromosome X 1 (Protein ALEX1) ARMCX1 ALEX1
Observed MW:	49-55 kDa
Uniprot ID :	Q9P291

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



mouse ovary tissue were subjected to SDS PAGE followed by western blot with FNab00588(ARMCX1 Antibody) at dilution of 1:300