

ARMCX2 antibody

Product Information

Catalog No.:	FNab00589
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	armadillo repeat containing, X-linked 2
Synonyms:	Armadillo repeat-containing X-linked protein 2 ARM protein lost in epithelial cancers on chromosome X 2 (Protein ALEX2) ARMCX2 ALEX2 KIAA0512
Observed MW:	66 kDa
Uniprot ID :	Q7L311

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:1000-1:4000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00589(ARM CX2 antibody) at dilution of 1:1000