

## CDNF antibody

### Product Information

Catalog No.:	FNab00593
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	arginine-rich, mutated in early stage tumors-like 1
Synonyms:	Cerebral dopamine neurotrophic factor ARMET-like protein 1 Conserved dopamine neurotrophic factor CDNF ARMETL1
Observed MW:	21 kDa
Uniprot ID :	Q49AH0

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:1000 ; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab00593( ARMETL1 Antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab00593( ARMETL1 Antibody) at dilution of 1:300