

ARNT antibody

Product Information

Catalog No.:	FNab00594
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	aryl hydrocarbon receptor nuclear translocator
Synonyms:	Aryl hydrocarbon receptor nuclear translocator (ARNT protein) Class E basic helix-loop-helix protein 2 (bHLHe2) Dioxin receptor, nuclear translocator Hypoxia-inducible factor 1-beta (HIF-1-beta, HIF1-beta) ARNT BHLHE2
Observed MW:	86 kDa
Uniprot ID :	P27540

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:2000
Image:	



IP Result of anti-ARNT, HIF1B (IP: FNab00594, 4ug; Detection: FNab00594 1:1000) with HepG2 cells lysate 3600ug.



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab00594 (ARNT, HIF1B antibody) at dilution of 1:1000