

ARSG antibody

Product Information

Catalog No.:	FNab00612
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Displays arylsulfatase activity at acidic pH with pseudosubstrates, such as p-nitrocatechol sulfate and also, but with lower activity, p-nitrophenyl sulfate and 4-methylumbelliferyl sulfate.

Immunogen information

Immunogen:	arylsulfatase G
Synonyms:	Arylsulfatase G (ASG) N-sulfoglucosamine-3-sulfatase ARSG KIAA1001
Observed MW:	62 kDa
Uniprot ID :	Q96EG1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab00612(ARSG antibody) at dilution of 1:800