

## CLMP antibody

### Product Information

Catalog No.:	FNab00619
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	adipocyte-specific adhesion molecule
Synonyms:	CXADR-like membrane protein Adipocyte adhesion molecule Coxsackie- and adenovirus receptor-like membrane protein (CAR-like membrane protein) CLMP ACAM ASAM
Observed MW:	45 kDa
Uniprot ID :	Q9H6B4

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:1000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00619( ASAM antibody at dilution of 1:600