

## ASB5 antibody

### Product Information

Catalog No.:	FNab00623
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May be a substrate-recognition component of a SCF-like ECS(Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. May play a role in the initiation of arteriogenesis(By similarity).

### Immunogen information

Immunogen:	ankyrin repeat and SOCS box-containing 5
Synonyms:	Ankyrin repeat and SOCS box protein 5 (ASB-5) ASB5
Observed MW:	36-40 kDa
Uniprot ID :	Q8WWX0

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab00623( ASB5 Antibody) at dilution of 1:300