

ASB8 antibody

Product Information

Catalog No.:	FNab00624
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be a substrate-recognition component of a SCF-like ECS(Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

Immunogen information


Immunogen:	ankyrin repeat and SOCS box-containing 8
Synonyms:	Ankyrin repeat and SOCS box protein 8 (ASB-8) ASB8
Observed MW:	32 kDa
Uniprot ID :	Q9H765

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab00624(ASB8 antibody) at dilution of 1:100



human heart tissue were subjected to SDS PAGE followed by western blot with FNab00624(ASB8 antibody) at dilution of 1:100