

## ASB9 antibody

### Product Information

Catalog No.:	FNab00625
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Substrate-recognition component of a SCF-like ECS(Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. Recognizes at least two forms of creatine kinase, CKB and CKMT1A.

### Immunogen information

Immunogen:	ankyrin repeat and SOCS box-containing 9
Synonyms:	Ankyrin repeat and SOCS box protein 9 (ASB-9) ASB9
Observed MW:	28-32 kDa
Uniprot ID :	Q96DX5

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab00625(ASB9 antibody) at dilution of 1:300