

ASF1B antibody

Product Information

Catalog No.:	FNab00634
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Histone chaperone that facilitates histone deposition and histone exchange and removal during nucleosome assembly and disassembly. Cooperates with chromatin assembly factor 1(CAF-1) to promote replication-dependent chromatin assembly. Does not participate in replication-independent nucleosome deposition which is mediated by ASF1A and HIRA. Required for spermatogenesis.

Immunogen information

Immunogen:	ASF1 anti-silencing function 1 homolog B
Synonyms:	Histone chaperone ASF1B Anti-silencing function protein 1 homolog B (hAsf1, hAsf1b) CCG1-interacting factor A-II (CIA-II, hCIA-II) ASF1B
Observed MW:	19 kDa
Uniprot ID :	Q9NVP2

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00634(ASF1B antibody) at dilution of 1:1000