

ATG2A antibody

Product Information

Catalog No.:	FNab00671
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for both autophagosome formation and regulation of lipid droplet morphology and dispersion.

Immunogen information


Immunogen:	ATG2 autophagy related 2 homolog A
Synonyms:	Autophagy-related protein 2 homolog A ATG2A KIAA0404
Observed MW:	215 kDa
Uniprot ID :	Q2TAZ0

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab00671(ATG2A Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00671(ATG2A Antibody) at dilution of 1:300