

## ACLY antibody

### Product Information

Catalog No.:	FNab00687
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	ATP citrate lyase
Synonyms:	ATP-citrate synthase ATP-citrate (pro-S-)-lyase (ACL) Citrate cleavage enzyme ACLY
Observed MW:	110-120 kDa
Uniprot ID :	P53396

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC, IP
Recommended dilution:	WB: 1:1000-1:4000; IP: 1:500-1:2000 for WB; IF: 1:10-1:100
Image:	



IP Result of anti-ATP citrate lyase (IP:FNab00687, 4ug; Detection:FNab00687 1:1000) with K-562 cells lysate 3200ug.



LO2 cells were subjected to SDS PAGE followed by western blot with FNab00687(ATP citrate lyase Antibody) at dilution of 1:1000